




Grandview

NonCycler Cow Program - Blanket Synch Group

CiDR 56hr CoSynch

	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
Nov 2011	7	8	9	10	11	12	13 CiDR IN + 1 mL GnRH at 8:00 a.m.-Blanket Synch Group - 
Nov 2011	14	15	16	17	18	19	20 CiDR OUT & 2 mL PG at 5:00 p.m.-Blanket Synch Group - 
Nov 2011	21 Breed on Heat + 1ml GnRH 	22 Breed on Heat + 1ml GnRH 	23 1 mL GnRH & Timed Breed at 9:00 a.m.-Blanket Synch Group - 	24	25	26	27
Nov 2011	28	29	30	1	2	3	4



GnRH - Give 1ml deep into the muscle using at least a 1" Needle 16ga or smaller with ratchet style injection gun



PG - Give 2ml deep into the muscle using at least a 1" Needle 16ga or smaller with ratchet style injection gun. PG should NOT be handled by pregnant women or asthmatics.



PG - Give 2ml deep into the muscle using at least a 1" Needle 16ga or smaller with ratchet style injection gun. PG should NOT be handled by pregnant women or asthmatics.



PMSG - Give 2ml deep into the muscle using at least a 1" Needle 16ga or smaller with ratchet style injection gun



GnRH - Give 2ml deep into the muscle using at least a 1" Needle 16ga or smaller with ratchet style injection gun



Ratchet Style Repeater Syringe
Adjustable from 1ml to 5ml